

INSTRUCTION FOR USE

Eosin-Floxin Alcoholic Solution

Catalog number	Size
05-B10020	1 x 150 mL
05-M10020	1 x 500 mL
05-10020/L	1 x 1 L

Expected aim

Product for the preparation of cyto-histological samples for optical microscopy.

Method :

The cytoplasmic staining Eosin, in combination with the nuclear staining Hematoxylin, is the most used method in the histopathological routine.

- 1- Section water
- 2- Mayer's hemalum 5 minutes
- 3- Tap water (or scott water for 1 minutes) 5 minutes
- 4- Washing in distilled water
- 5- Ethanol 80% 1 minute
(optional step to preserve the alcoholic content of the solution of eosin-phloxine)
- 6- Eosin-floxin alcoholic Solution 2-3 minutes
- 7- Ethanol 95% x 2 1 minute
- 8- Ethanol 95% 1 minute
- 9- Absolute ethanol 2 minutes
- 10- Absolute ethanol 2 minutes
- 11- Clearing agent and mount

Warning and Precaution

The product must be used exclusively by specialized technical operators.

Carefully read the information on the classification of dangerous substances on the label. Always refer to the safety data sheet where are available the information on the risk presented by the mixture, the precautionary measures during use, the measures first aid and the event of accidental release.

Do not use if the primary container is damaged.

Contraindication

Obtain special instructions before use.

Wear protective gloves / clothing and eye / face protection.

Wash hands thoroughly after handling

IF INHALED: remove victim to fresh air and keep at rest position comfortable for breathing.

Call a POISON CENTER or doctor/physician if you feel unwell.