

INSTRUCTION FOR USE

Product name: Bio Cassettes II

Code	Description
07-8100	Embedding cassettes, white - 2000 pcs
07-8110	Embedding cassettes, orange - 2000 pcs
07-8120	Embedding cassettes, blue - 2000 pcs
07-8130	Embedding cassettes, yellow - 2000 pcs
07-8140	Embedding cassettes, lilac - 2000 pcs
07-8150	Embedding cassettes, pink - 2000 pcs
07-8160	Embedding cassettes, green - 2000 pcs

Expected aim

Product for the preparation of cyto-histological samples for optical microscopy.

Purpose of use

To contain histology specimen in the tissue processor

Mode of use

- Put the specimen inside the cassette;
- Close the cassette.