

## INSTRUCTION FOR USE

### Dehyol Absolute

Catalog number	Size
06-10077E	1 x 2,5 L
06-10077Q	4 x 2,5 L
06-10077F	1 x 5L

### Expected aim

Product for the preparation of cyto-histological samples for optical microscopy.

### Method :

Put the solvent into staining dishes or automatic stainer / tissue processors according to staining or processing protocol.

### Warning and Precaution

The product must be used exclusively by specialized technical operators.  
Carefully read the information on the classification of dangerous substances on the label. Always refer to the safety data sheet where are available the information on the risk presented by the mixture, the precautionary measures during use, the measures first aid and the event of accidental release.  
Do not use if the primary container is damaged.

### Contraindication

Obtain special instructions before use.  
Wear protective gloves / clothing and eye / face protection.  
Wash hands thoroughly after handling.  
Call a POISON CENTER or doctor/physician if you feel unwell.  
IF INHALED: remove victim to fresh air and keep a rest in position comfortable for breathing.