

INSTRUCTION FOR USE

Oil Red O

| <u>Catalog number</u> | <u>Size</u> |
|-----------------------|-------------|
| 04-220923 | 100 test |

Expected aim

Product for the preparation of cyto-histological samples for optical microscopy.

Method :

- 1) Bring the sections to the distilled water.
- 2) Put the reagent A in a Coplin Jar and immerse the slide for 20 minutes.
- 3) Wash quickly in tap water.
- 4) Drain the slide and put on the section 10 drops of reagent B: leave to act 30 seconds.
- 5) Wash in running tap water for 3 minutes.
- 6) Drain and mount in aqueous medium.

Warning and Precaution

The product must be used exclusively by specialized technical operators.
Carefully read the information on the classification of dangerous substances on the label. Always refer to the safety data sheet where are available the information on the risk presented by the mixture, the precautionary measures during use, the measures first aid and the event of accidental release.

Do not use if the primary container is damaged.

Contraindication

Obtain special instructions before use.

Wear protective gloves / clothing and eye / face protection.

Keep away from heat/sparks/open flames/hot surface. No Smoking.

IF ON SKIN (on hair): remove/take off immediately all contaminated clothing. Rinse skin with water/shower.