

BIOCASSETTE II



CODE	COLOUR	PACKAGING
07-8100	White	2000 pcs (3 dispensers: 2 with 1000 cassettes and 1 with 2000 lids)
07-8110	Orange	2000 pcs (3 dispensers: 2 with 1000 cassettes and 1 with 2000 lids)
07-8120	Blue	2000 pcs (3 dispensers: 2 with 1000 cassettes and 1 with 2000 lids)
07-8130	Yellow	2000 pcs (3 dispensers: 2 with 1000 cassettes and 1 with 2000 lids)
07-8140	Lilac	2000 pcs (3 dispensers: 2 with 1000 cassettes and 1 with 2000 lids)
07-8150	Pink	2000 pcs (3 dispensers: 2 with 1000 cassettes and 1 with 2000 lids)
07-8160	Green	2000 pcs (3 dispensers: 2 with 1000 cassettes and 1 with 2000 lids)



In Vitro Diagnostic – medical device



Manufacturer: Bio-Optica Milano S.p.A.



Disposable product

Date of emission: 25/03/2020

Rev. 03

Standard cassettes for histological samples embedding with separated lid.

TECHNICAL SPECIFICATIONS

Cassettes for the embedding of histological samples with separated lid, with a double closure system to guarantee a perfect lock-up during the processing passage. The sample's identification is possible thanks to the wizen surface on three cassette's side. These surfaces are fitted to the writing with pencil. For a permanent marking, we suggest to use Histology pen (code 11-50 or 11-400). These cassettes are compatible with Bio-Optica filing system (Istobloc, Colorteca and big archives).

Biocassettes are manufactured with our best technical skills, therefore, due to the complexity of the details themselves and the type of material used, is to be considered a physiological quantity of defective cassettes equal to 5,000 PPM (Parts Per Million), corresponding to a percentage of 0.5%. Within this amount, Bio-Optica does not accept complaints and / or requests for warranty replacement.

Technical data	Chemical composition	Acetalic polymer (polyoximethylene polyacetal). The acetalic polymer is stable to high temperatures, microwave operations, solvents, acids and alkalis; therefore this material is optimal for the processing and the embedding passages in the histological process.
	Internal dimensions	30 x 25 x 5 mm
	External dimensions	40 x 28 x 6 mm
	Holes dimensions	(5.8 ± 0.2) x (1 ± 0.1) mm
	Surface angle	45°
Packaging	Primary container	Carton dispenser with frontal opening to ease the extraction of cassettes
	Secondary container	Carton box
Storage	Storing conditions	None
	Validity	Not applicable for this kind of product
Warning and precautions	Product classification	For specialized technical operators. Product not classified as dangerous.
	Disposal	Dispose in conformity with local laws.

Date of emission: 25/03/2020

Rev. 03